

HT-1
Pencil Couter

HT-10
Pencil Couter

HT-5
Pencil Couter

ESU Pencil
HUATONG Brand

Disposable Electro Surgical
Pencil Couter


The HUATONG Single Use Electrosurgical pencil. The electrode is fitted with a hexagonal connection, securing the electrode in the required position, it can be used with any electrosurgical generator available.

- Comfortable design and light weight
- Colour market buttons
- Hexagonal electrode connection
- 3 meter cabl, 2,36mm diameter electrode
- 3 pin US-connector

Holster is used for place ESU pencil, protect surgery safety. The tip cleaner (tip cleaner sponge) used for cleaning the tissue left on the eletrosurgical pencil tip, sticky backing for fixing on the surgical table facilitates cleaning the tip.

Details:

Length of handle	170mm (Only for HP05) 156MM (The others)	Package of each piece: Packed with on PE plastic bag + one medical paper pouch	
Length of Blade	70mm (40mm outside of pencil)		
Dimaeter of Blade	2.36mm	PACKAGE OF EACH CARTON	
Length of wire	10ft (3M)	QTY	100pcs/box
Maximal output current	1A	Size of box	38 x 30 x 23 cm
Frenquency confine	0-1.0 MHZ	Gross weight	7.7 kg
Current resistance	200 Ω	Net weight	7.0 kg
Sterilized method	Strrile by EO	20FT Container can fit	100.000 PCS
Term of sterility	3 years	40FT Container can fit	240.000 PCS

HT-7 ESU Pencil**HT-9 ESU Pencil****HT-5 ESU Pencil****Spesification**

- ☑ 153mm-ABS pencil appearance, foot control
- ☑ 7cm (40mm exposed length) Stainless Steel tip
- ☑ 2.36mm-diameter of the tip
- ☑ 3m-PVC medical cable, ϕ 6.3mm one pin jack phone plug
- ☑ Maximal output current: 1A
- ☑ Frequency confine: 0-1.0MHZ
- ☑ Sterile by EO.

Spesification

- ☑ 153mm-ABS pencil appearance, smoke-evacuation
- ☑ 4-3m PVC medical cable
- ☑ Maximal output current: 1A
- ☑ Frequency confine: 0-1.0MHZ
- ☑ Sterile by EO
- ☑ Current resistance: 200 Ω
- ☑ Term of sterility: 3 years
- ☑ Package of each piece: Package with 1PE plastic bag+1 medical paper pouch.

Disposable Neutral Electrode**HT-PA1-WV and HT-PA2-WV****HT-PA1T and HT-PA2T****HT-PA1T-WV and HT-PA2T-WV****HT-PA1 and HT-PA2****HT-PB1 | HT-PB2 | HT-PB2-WV****Spesification**

- ☑ Suitable for different kinds of surgery
- ☑ Special adhesive layer proof firm attachment, safe and comfortable
- ☑ Provide adult, pediatric, and infants series
- ☑ Monopolar type and Bipolar type
- ☑ Pre-wired and non-wired both OK
- ☑ Disposable Neutral Electrode with Cable (Grounding Pad/Dispersive Electrode/Neutral Pad).

Disposable Electrosurgical Pad, also named ESU Pad, Grounding Pad, Neutral Plate, Patient Plate, and Dispersive Electrode etc., is approved by CE from TUV. We are not only offers basic Electrosurgical Pads of mono-polar and bipolar (solid and split foil) in upright and horizontal design fit for adult.



PT. GOLDEN GLOBE MEDICA
providing a better and reliable solution

**Head Office**

Sentra Industri Terpadu Blok E2 No. 25, Pantai Indah Kapuk, Jakarta Utara 14470
Jakarta - Indonesia

Telp : +62 21 5698 2721 | +62 21 5698 3275 Telp : +62 21 350 3011
Fax : +62 21 5698 2720 | +62 21 5698 3276 Fax : +62 21 350 2667

Representative Office

● Yogyakarta
Sangrahan No. 037 RT/RW.08/09 Banyuraden Gamping, Sleman
Yogyakarta 55293
Telp/ Fax : +62 274 530 4195